

EarlyCDT[®]—Lung



A simple blood test to aid in indeterminate pulmonary nodule risk assessment.



Complement to CT scans



Measures a panel of autoantibodies



Characterizes risk of indeterminate nodules




May reclassify nodule to intervention risk



Earlier intervention and better outcomes

EarlyCDT—Lung Reclassifies Malignancy Risk in Nodules with 10%–65% Initial Risk



A nodule risk calculator combining the Swensen/Mayo model with **EarlyCDT—Lung** is available at <http://oncimmune.com/nodule-calculator>
 Risk calculations: Healey GF, Macdonald IK, Reynolds C, et al. Tumor-Associated Autoantibodies: Re-Optimization of **EarlyCDT—Lung** Diagnostic Performance and its Application to Indeterminate Pulmonary Nodules. *J Cancer Ther.* 2017; 8:506–517.



Strong clinical data demonstrating high specificity and high PPV (Rule-in Test)



Highly cost-effective and may be billed to Medicare and commercial insurance companies



Test reports show complete information at a glance

Nodule Risk Assessment Patient Report

Note: This test is not intended for patients who have a previous diagnosis of cancer (except basal cell carcinoma).

Patient	Specimen Information	Clinician
PATIENT 5	T	Dr. H.L.7 Physician
Sex: F	C	Dr. H.L.7 Clinic
D.O.B: 7/4/1962	F	Sesame St.
Age: 55	Report date: 8/31/2017	De Soto, KS 66018

EarlyCDT—Lung Test Result: HIGH LEVEL

Test	Result (RU)	No Significant Level of Autoantibodies Detected	Moderate Level	High Level
CAGE autoantibody	<2.76	-X-	---	---
GBU4-5 autoantibody	1.06	-X-	---	---
NY-ESO-1 autoantibody	<1.01	-X-	---	---
p53 autoantibody	<3.09	-X-	---	---
SOX-2 autoantibody	5.62	---	---	-X-
MAGE A4 autoantibody	<3.91	-X-	---	---
HuD autoantibody	<3.99	-X-	---	---

Clinical Utility

The ACCP guidelines^a recommend assessing the risk of malignancy of a pulmonary nodule, e.g., with the Swensen/Mayo nodule malignancy risk calculator^b, available at oncimmune.com/nodule-calculator. The calculated risk can be divided into three categories and the patient managed accordingly.

EarlyCDT—Lung facilitates further risk characterization to assist with triaging difficult to assess nodules.^{c,d}

<5% risk of lung cancer* VERY LOW RISK	High or Moderate EarlyCDT—Lung test result: risk raised from very low risk to low to moderate risk.
5%–65% risk of lung cancer* LOW to MODERATE RISK	High EarlyCDT—Lung test result: risk raised to high risk if pre-test risk >10%. Moderate EarlyCDT—Lung test result: risk raised to high risk if pre-test risk >45%; otherwise, consider patient at increased moderate risk.
>65% risk of lung cancer* HIGH RISK	Occasional use of EarlyCDT—Lung following biopsy or bronchoscopy where further risk evaluation is deemed of value.

* Risk categories according to the ACCP guidelines.¹

Interpretive Comments

A **High Level** result is reported when any one or more autoantibodies in the **EarlyCDT—Lung** panel are above the high cut-off value. For a nodule with a pre-test risk of >10%, a High Level **EarlyCDT—Lung** result will move the nodule to high risk (>65%). Consider changing the patient's treatment pathway to that recommended by guidelines for a nodule at high risk of malignancy.

References

a) Gould MK, et al. *Chest* 2013; 143(5):e93S–e120S.
 b) Swensen SJ, et al. *Arch Intern Med*. 1997; 157:849–853.
 c) Massion pp, et al. *J Thorac Oncol*. 2017; 12(3):578–583.
 d) This test was developed and its performance characterized by a CLIA accredited, high complexity laboratory and is performed by Joseph P. McConnell, PhD, DABCC, FACB, Clinical Laboratory Director at the University of Kansas Medical Center.

www.oncimmune.com

Questions? Call 1-888-583-9030

What do these biomarker levels have to do with liver cancer?

Some individuals with liver cancer have been found to have elevated levels of one or more of these biomarkers.^{a–f} Autoantibodies have been shown to be present in the blood up to four years prior to cancer diagnosis by imaging.^{g,i} Early detection of liver cancer has been shown to increase the potential for an improved outcome.^h

References

a) Middleton CH, et al. *PLoS ONE* 2014; 9(8):e103867.
 b) Hong Y, et al. *World J Hepatol*. 2015; 7(11):1581–1585.
 c) Chien Y, et al. *Cancer Letters* 2010; 289:32–39.
 d) Liu H, et al. *Cancer Epidemiol*. 2012; 26(1):82–88.
 e) Tatarinov IS. *Vopr Med Khim*. 1964; 10:90–91.
 f) Hu B, et al. *Int J Mol Sci*. 2013; 14(12):23559–23580.
 g) Trivers GE, et al. *Clin Cancer Res*. 1996; 2:1767–1775.
 h) Li Y, et al. *Int J Cancer*. 2005; 114:157–160.
 i) Zhong L, et al. *J Thorac Oncol*. 2006; 1:513–519.
 j) Jett J, et al. *J Thorac Oncol*. 2017; 12(11):S2170.
 k) Zhang BH, et al. *J Cancer Res Clin Oncol*. 2004; 130(7):417–422.

www.oncimmune.com

Questions? Call 1-888-583-9030

Individual antibody results

Detailed breakdown of cutoff levels and further education on results for patients

Link to risk calculator

Cut-off value	Moderate Level result	High Cut-off value	High Level result
		8.696	
6.865	9.696	6.574	6.774
7.972	8.172	6.286	6.606
6.515	6.715	7.826	8.026
00 ng/mL	221 ng/mL		

Standard cut-off applied when AFP is used as a single biomarker test to monitor for recurrence after liver transplant.