



How to use *EarlyCDT—Lung*



* these may vary according to geographic region

How does *EarlyCDT—Lung* work?

EarlyCDT—Lung measures blood levels of a panel of seven autoantibodies to tumour-associated antigens that are linked to lung cancer. The seven autoantibodies have been shown to be elevated for all types of lung cancer, and from the earliest stage of the disease.^{1,2} Unlike the tumour antigens themselves, the autoantibody levels can be measured easily and accurately, based upon the signal magnification created by the body's immune response to cancer. The test runs on a simple enzyme-linked immunosorbent assay (ELISA) platform, which is widely available in hospital laboratories around the world and in Oncimmune's CLIA laboratory in Kansas, US.



More than **155,000 tests** have already been performed for patients worldwide, which represents tests ordered by more than **2,000 clinicians**. *EarlyCDT—Lung* detects all types and stages of lung cancer and has led to the detection of numerous early stage lung cancers. *EarlyCDT—Lung* is now available in many countries. Please see the list of test providers at <http://oncimmune.com/distributors/>

- Noone AM, Howlader N, Krapcho M, et al. *SEER Cancer Statistics Review, 1975–2015*, National Cancer Institute, Bethesda, MD, https://seer.cancer.gov/csr/1975_2015/, based on November 2017 SEER data submission, posted to the SEER web site, April 2018.
- U.S. Preventive Services Task Force. Final Recommendation Statement: Lung Cancer: Screening. December 2016. <http://uspreventiveservicestaskforce.org/Page/Document/UpdateSummaryFinal/lung-cancer-screening>.
- The National Lung Screening Trial Research Team. Reduced lung-cancer mortality with low-dose computed tomographic screening. *N Engl J Med*. 2011; 365:395–409.
- Jett J, Healey G, MacDonald I, et al. Determination of the detection lead time for autoantibody biomarkers in early stage lung cancer using the UKCTOCS cohort. *J Thorac Oncol*. 2017;12(11):S2170.
- Zhong L, Coe SP, Stromberg AJ, et al. Profiling tumor-associated antibodies for early detection of non-small cell lung cancer. *J Thorac Oncol*. 2006;1(6):513–519.
- Sullivan F, Dorward A, Mair F, et al. P2.06-038. An RCT of the detection of autoantibodies to tumour antigens in lung cancer using the EarlyCDT-Lung test in Scotland (ECLS) in 12,208 study subjects. *J Thorac Oncol*. 2017;12(1):S1095.
- Lung Cancer Alliance, Types of Lung Cancer, <http://lungcanceralliance.org/get-information/types-of-lung-cancer.html>.
- Chapman CJ, Healey GF, Murray A, et al. EarlyCDT®-Lung test: improved clinical utility through additional autoantibody assays. *Tumor Biol*. 2012;33(5):1319–26.
- Boyle P, Chapman CJ, Holdenrieder S, et al. Clinical validation of an autoantibody test for lung cancer. *Ann Oncol*. 2011;22(2):383–389.
- Healey GF, Lam S, Boyle P, et al. Signal stratification of autoantibody levels in serum samples and its application to the early detection of lung cancer. *J Thorac Dis*. 2013;5(5):618–625.
- National Cancer Institute. Breast Cancer Surveillance Consortium: Evaluating Screening Performance in Practice. *NIH Publication No. 04-5490*. Bethesda, MD: National Cancer Institute, National Institutes of Health, U.S. Department of Health and Human Services, April 2004.
- Imperiale TF, Ransohoff DF, Itzkowitz SH et al. Multitarget stool DNA testing for colorectal-cancer screening. *N Engl J Med*. 2014; 370(14):1287–1297
- Positive Predictive Value – the number of positive test results required to detect a cancer.
- Exception: basal cell carcinoma. See *EarlyCDT—Lung* FAQs

+44 (0)115 823 1869

www.oncimmune.com

Oncimmune®



EarlyCDT—Lung

A simple blood test to aid in the risk assessment and early detection of lung cancer in high-risk patients.

Tel +44 (0) 115 8231869
Fax +44 (0) 115 8231958
earlycdt@oncimmune.co.uk
www.oncimmune.com

Oncimmune Limited
Clinical Sciences Building,
City Hospital, Hucknall Road,
Nottingham NG5 1PB